EnVision™ FLEX Mini Kit, High pH, (Dako Autostainer/Autostainer Plus)

Code K8024

3rd edition/ 3^{ème} édition/ 3. Ausgabe

The kit contains reagents sufficient for 125-190 tests. For Dako Autostainer/Autostainer Plus.

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Page/ Page/ Seite

ENGLISH

Intended Use	4
Summary and Explanation	4
Reagents	5
A. Materials provided	5
B. Materials required but not provided	5
C. Optional reagents	6
Precautions	6
Storage	8
Reagent Preparation	8
Specimen Collection and Preparation	9
Procedure	9
A. Pre-treatment procedure	9
B. Staining procedure	10
Quality Control	11
Interpretation of Results	11
Limitations	11
Explanation of symbols	29

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Intended Use

For in vitro diagnostic use.

Dako EnVision™ FLEX detection system is intended for use in immunohistochemistry together with Dako Autostainer/Autostainer Plus instruments. The system detects primary mouse and rabbit antibodies and the reaction is visualized by EnVision™ FLEX DAB+ Chromogen. If used with EnVision™ FLEX+ Rabbit (LINKER) (Code K8019) or EnVision™ FLEX+ Mouse (LINKER) (Code K8022) signal amplification of primary rabbit antibodies or primary mouse antibodies, respectively, can be achieved. The EnVision™ FLEX reagents are intended for use on formalin-fixed, paraffin-embedded tissue sections.

Summary and Explanation

Dako EnVision™ FLEX detection system, has been designed to be a flexible system and gives an optimal staining on Dako Autostainer/Autostainer Plus instruments, when using the protocol recommended in this package insert.

Prior to staining formalin-fixed, paraffin-embedded tissue sections should be subjected to deparaffinization and hydration followed by heat-induced epitope retrieval (HIER) using the target retrieval method specified in the package insert for the primary antibody. Please also see Procedure section for 3-in-1 specimen preparation. Code K8024 includes EnVision™ FLEX Target Retrieval Solution, High pH (50x) (DM828). Alternatively, EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8004) or EnVision™ FLEX Target Retrieval Solution, Low pH (50x) (Code K8005) is recommended. Some primary antibodies require enzymatic pre-treatment of tissue for optimal staining instead of HIER.

Endogenous peroxidase should be blocked with EnVision™ FLEX Peroxidase-Blocking Reagent (DM821) included in the kit. Due to an effective washing procedure and the presence of stabilizing proteins in the Dako reagents, extra blocking steps to reduce non-specific background staining are unnecessary.

EnVision™ FLEX Wash Buffer (20x) (DM831) included in the kit or available as optional reagent (Code K8007) is recommended.

Primary antibodies are not provided with the kit. We recommend the use of Dako FLEX Ready-to-Use Primary Antibodies or Dako concentrated Primary Antibodies. EnVision™ FLEX Antibody Diluent (Code K8006) is recommended for the dilution of Dako concentrated Primary Antibodies.

EnVision™ FLEX+ Rabbit (LINKER) (Code K8019) may be applied for an optional signal amplification of rabbit primary antibodies. For signal amplification when using mouse primary antibodies, EnVision™ FLEX+ Mouse (LINKER) (Code K8022 or included in Code K8012) may be applied. As a guideline: EnVision™ FLEX+ Rabbit (LINKER) provides 2-3 fold signal amplification and EnVision™ FLEX+ Mouse (LINKER) provides 4-5 fold signal amplification.

Dako EnVision™ FLEX /HRP detection reagent (DM822) in the kit consists of a dextran backbone to which a large number of peroxidase (HRP) molecules and secondary antibody molecules have been coupled. A unique chemistry is used for the coupling reaction, which permits the binding of up to 100 HRP molecules and up to 20 antibody molecules per backbone.

The substrate system in the kit consists of two components: EnVision™ FLEX DAB+ Chromogen (DM827), a concentrated diaminobenzidine (DAB) solution, and EnVision™ FLEX Substrate Buffer (DM823) containing hydrogen peroxide. Before use EnVision™ FLEX DAB+ Chromogen must be diluted in EnVision™ FLEX Substrate Buffer. The substrate system produces a crisp brown end product at the site of the target antigen.

EnVision™ FLEX Hematoxylin (Code K8018) is recommended for counterstaining. The reagent provides a clear blue, nuclear staining.

The stained tissue sections may be mounted with either aqueous or organic-solvent-based mounting medium.

Reagents

A. Materials provided

DM821

EnVision™ FLEX
PEROXIDASE-BLOCKING
REAGENT

EnVision™ FLEX Peroxidase-Blocking Reagent

3 x 13 mL, ready-to-use

Phosphate buffer containing hydrogen peroxide, 15 mmol/L NaN₃ and detergent.

DM822

EnVision™ FLEX
/HRP

EnVision™ FLEX /HRP

3 x 13 mL, ready-to-use

Dextran coupled with peroxidase molecules and goat secondary antibody molecules against rabbit and mouse immunoglobulins. In buffered solution containing stabilizing protein and preservative.

DM827

EnVision™ FLEX
DAB+ CHROMOGEN

EnVision™ FLEX DAB+ Chromogen

1 x 3 mL

3,3'-diaminobenzidine tetrahydrochloride in organic solvent.

The color of this reagent may vary from strong violet to colorless without having any influence on the performance of the kit.

DM823

EnVision™ FLEX
SUBSTRATE BUFFER

EnVision™ FLEX Substrate Buffer

7 x 13 mL

Buffered solution containing hydrogen peroxide and preservative.

DM828

EnVision™ FLEX
TARGET RETRIEVAL SOLUTION
HIGH pH (50x)

EnVision™ FLEX Target Retrieval Solution, High pH (50x)

3 x 30 mL, 50x concentrated

Tris/EDTA buffer, pH 9

DM831

EnVision™ FLEX
WASH BUFFER (20x)

EnVision™ FLEX Wash Buffer (20x)

2 x 1 L, 20x concentrated

Tris-buffered saline solution containing Tween 20, pH 7.6 (±0.1).

B. Materials required but not provided

Dako Autostainer/Autostainer Plus instrument

Dako PT Link

Dako Proteinase K, Ready-to-Use (if necessary)

Dako FLEX Ready-to-Use Primary Antibodies or suitably diluted concentrated primary rabbit or mouse antibodies from Dako

Dako FLEX Ready-to-Use Rabbit or Mouse Universal Negative Control or suitable negative control reagent for the primary antibody

Dako instrumentation utensils

Microscope slides, e.g. FLEX IHC Microscope Slides

General laboratory reagents for deparaffinization of paraffin-embedded tissue sections
Mounting medium (aqueous or organic-solvent-based) and coverslips

C. Optional reagents

K8004
DM828

EnVision™ FLEX
TARGET RETRIEVAL SOLUTION
HIGH pH (50x)

EnVision™ FLEX Target Retrieval Solution, High pH (50x)
3 x 30 mL, 50x concentrated
Tris/EDTA buffer, pH 9

K8005
DM829

EnVision™ FLEX
TARGET RETRIEVAL SOLUTION
LOW pH (50x)

EnVision™ FLEX Target Retrieval Solution, Low pH (50x)
3 x 30 mL, 50x concentrated
Citrate buffer, pH 6.1

K8006
DM830

EnVision™ FLEX
ANTIBODY DILUENT

EnVision™ FLEX Antibody Diluent
120 mL, ready-to-use
Tris buffer, pH 7.2, containing 15 mmol/L NaN₃, and protein.

K8007
DM831

EnVision™ FLEX
WASH BUFFER (20x)

EnVision™ FLEX Wash Buffer (20x)
1 L, 20x concentrated
Tris-buffered saline solution containing Tween 20, pH 7.6 (±0.1).

K8018
DM826

EnVision™ FLEX
HEMATOXYLIN

EnVision™ FLEX Hematoxylin
10 x 13 mL, ready-to-use
Aqueous solution of hematoxylin.

K8019
DM825

EnVision™ FLEX+
RABBIT (LINKER)

EnVision™ FLEX+ Rabbit (LINKER)
3 x 13 mL, ready-to-use
Buffered solution containing stabilizing protein and an antimicrobial agent.
Contains reagent sufficient for 110-165 tests based on the use of 200-300 µL per test.

K8022
SM824

EnVision™ FLEX+
MOUSE (LINKER)

EnVision™ FLEX+ Mouse (LINKER)
3 x 13 mL, ready-to-use
Buffered solution containing stabilizing protein and an antimicrobial agent.
Contains reagent sufficient for 130-200 tests based on the use of 200-300 µL per test.

Precautions

1. For professional users.
2. EnVision™ FLEX /HRP, EnVision™ FLEX+ Mouse (LINKER) and EnVision™ FLEX+ Rabbit (LINKER) contain material of animal origin and it cannot be excluded that trace amounts of human material could be present due to manufacturing procedures. As with any product derived from biological sources, proper handling should be used.

3. Do not expose EnVision™ FLEX DAB+ Chromogen, EnVision™ FLEX Substrate Buffer or the EnVision™ FLEX Substrate Working Solution to strong light during the procedure.
4. EnVision™ FLEX Target Retrieval Solution, High pH (50x) and EnVision™ FLEX Target Retrieval Solution, Low pH (50x) contain 5-10% Nonoxinol and are labeled:

Warning

H319	Causes serious eye irritation.
H411	Toxic to aquatic life with long lasting effects.
P280	Wear eye or face protection.
P273	Avoid release to the environment.
P264	Wash hands thoroughly after handling.
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337 + P313	If eye irritation persists: Get medical attention.
P501	Dispose of contents and container in accordance with all local, regional, national and international regulations.

5. For EnVision™ FLEX Wash Buffer (20x) contains 10-30% 2-Amino-2-(hydroxymethyl)propane-1,3-diol hydrochloride and is labeled:

Warning

H319	Causes serious eye irritation.
H315	Causes skin irritation.
P280	Wear protective gloves. Wear eye or face protection.
P264	Wash hands thoroughly after handling.
P302 + P352 + P362-2 + P363	IF ON SKIN: Wash with plenty of soap and water. Take off contaminated clothing. Wash contaminated clothing before reuse.
P332 + P313	If skin irritation occurs: Get medical attention.
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337 + P313	If eye irritation persists: Get medical attention.

6. For EnVision™ FLEX Peroxidase-Blocking Reagent, a safety data sheet is available for professional users on request.
7. EnVision™ FLEX DAB+ Chromogen contains 1-5% 3,3'-diaminobenzidine tetrahydrochloride, and is labeled:

Danger

H350	May cause cancer.
H341	Suspected of causing genetic defects.
P201	Obtain special instructions before use.
P202	Do not handle until all safety precautions have been read and understood.
P280	Wear protective gloves. Wear eye or face protection. Wear protective clothing.
P308 + P313	IF exposed or concerned: Get medical attention.
P405	Store locked up.
P501	Dispose of contents and container in accordance with all local, regional, national and international regulations.

8. Optional reagent EnVision™ FLEX Antibody Diluent (Code K8016) contains sodium azide (NaN_3), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
9. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
10. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store EnVision™ FLEX and EnVision™ FLEX+ reagents at 2–8 °C.

Store EnVision™ FLEX Hematoxylin (Code K8018) at room temperature. Do not use after expiration date stamped on the vials. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. The prepared Substrate Working Solution should be stored in the dark at 2–8 °C and used within 5 days. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the product is suspected, contact Dako Technical Services.

Reagent Preparation

A.1 EnVision™ FLEX Target Retrieval Solution, High pH (50x) (DM828, Code K8004)

Dilute a sufficient quantity of EnVision™ FLEX Target Retrieval Solution, High pH (50x) 1:50 using distilled or deionised water for the staining procedure that is planned. Unused diluted solution may be stored at 2–8 °C for one month. If Dako PT Link is used for pre-treatment, dilution can be performed by emptying the content of the Target Retrieval Solution (50x) vial into the Dako PT Link tank and adding distilled or deionised water to the marked line (Fill Line). Discard solution if cloudy in appearance. When used in PT Link for 3-in-1 specimen preparation procedure, the diluted solution can be used 3 times within a 5 day period if stored at room temperature.

A.2 EnVision™ FLEX Target Retrieval Solution, Low pH (50x) (DM829, Code K8005)

Dilute a sufficient quantity of EnVision™ FLEX Target Retrieval Solution, Low pH (50x) 1:50 using distilled or deionised water for the staining procedure that is planned. If Dako PT Link is used for pre-treatment, dilution can be performed by emptying the content of the Target Retrieval Solution (50x) vial into the Dako PT Link tank and adding distilled or deionised water to the marked line (Fill Line). Unused diluted solution may be stored at 2–8 °C for one month. Discard solution if cloudy in appearance. When used in PT Link for 3-in-1 specimen preparation procedure, the diluted solution can be used 3 times within a 5 day period if stored at room temperature.

A.3 EnVision™ FLEX Wash Buffer (20x) (DM831, Code K8007)

Dilute a sufficient quantity of EnVision™ FLEX Wash Buffer (20x) 1:20 using distilled or deionised water for the staining procedure that is planned. Dilute by adding the concentrated Wash Buffer to a pre-measured quantity of distilled or deionised water to minimize foaming. Stir gently until the diluted solution appears homogenous. Unused diluted solution may be stored at 2–8 °C for one month. Discard solution if cloudy in appearance. When used in PT Link for 3-in-1 specimen preparation procedure, the diluted solution can be used 3 times within a 5 day period if stored at room temperature.

A.4 EnVision™ FLEX Substrate Working Solution

The DAB-containing EnVision™ FLEX Substrate Working Solution is prepared by mixing thoroughly 1 drop EnVision™ FLEX DAB+ Chromogen (DM827) per 1 mL EnVision™ FLEX Substrate Buffer (DM823). Use EnVision™ FLEX Substrate Working Solution within 5 days (store in the dark at 2–8 °C).

Specimen Collection and Preparation

The specimens may be formalin-fixed, paraffin-embedded tissue sections.

Fixation time is dependent on fixative and tissue type/thickness. For example, tissue blocks with a thickness of 3–4 mm should be fixed in neutral-buffered formalin for 18–24 hours.

The optimal thickness of paraffin-embedded sections is approximately 4 µm.

The specimens should be mounted on microscope slides, e.g. FLEX IHC Microscope Slides (Code K8020). The sections should be mounted on the slides as flat and wrinkle-free as possible. Too many wrinkles will have an impact on the staining results.

NOTE: The microscope slides must have a width suitable for the Dako Autostainer/Autostainer Plus instrument. Please refer to the Operator's Manual for the individual Dako instrument for definition of usable microscope slides.

Paraffin sections should be mounted from a pre-heated water bath containing distilled or deionized water. The water bath should contain no additives (such as gelatin, poly-L-lysine etc.). Sections should be dried by heating, generally at a temperature not above 60 °C for up to 60 minutes. To ensure proper adherence of sections to slides it is important to drain the water from beneath the sections prior to the oven drying process.

Procedure

Dako Autostainer/Autostainer Plus instruments use techniques based on different principles to obtain an optimal staining result. Before running protocols on your Dako Autostainer/Autostainer Plus instrument, please read carefully the **Operator's Manual for the dedicated Dako instrument**.

In order to ascertain reproducible epitope retrieval it is recommended always to load the slide holder fully with slides. This ensures an identical heating of sections in every run.

A. Pre-treatment procedure

Recommended 3-in-1 specimen preparation procedure using PT Link:

Deparaffinization, rehydration and heat-induced epitope retrieval (HIER) can be performed on formalin-fixed paraffin-embedded tissue sections using the 3-in-1 procedure:

1. Prepare a working solution by diluting the EnVision™ FLEX Target Retrieval Solution concentrate 1:50 in distilled or deionized water.
2. Fill PT Link tanks with sufficient quantity (1.5 L) of working solution to cover the tissue sections.
3. Set PT Link to pre-heat the solution to 65 °C.
4. Immerse the mounted, formalin-fixed, paraffin-embedded tissue sections into the pre-heated EnVision™ FLEX Target Retrieval Solution (working solution) in PT Link tanks and incubate for 20-40 minutes at 97 °C. The optimal incubation time should be determined by the user.
5. Leave the sections to cool in PT Link to 65 °C.
6. Remove each Autostainer slide rack with the slides from the PT Link tank and immediately dip slides into a jar/tank (e.g., PT Link Rinse Station, Code PT109) with diluted, room temperature EnVision™ FLEX Wash Buffer (20x).
7. Leave slides in the diluted, room temperature EnVision™ FLEX Wash Buffer (20x) for 1-5 minutes.
8. Place slides on a Dako Autostainer/Autostainer Plus instrument and proceed with staining. The sections should not dry out during the treatment or during the immunohistochemical staining procedure.
9. After staining it is recommended to perform dehydration, clearing and permanent mounting.

NOTE: When used in PT Link for 3-in-1 specimen preparation procedure, the diluted Target Retrieval Solution and diluted Wash Buffer can be used 3 times within a 5 day period if stored at room temperature.

HIER procedure using PT Link:

After deparaffinization and hydration to buffer (water), tissue sections should be subjected to heat-induced epitope retrieval (HIER):

1. Prepare a working solution by diluting the EnVision™ FLEX Target Retrieval Solution (50x) concentrate 1:50 in distilled or deionized water.
2. Fill PT Link tanks with sufficient quantity (1.5 L) of working solution to cover the tissue sections.
3. Set PT Link to pre-heat the solution to 65 °C.
4. Manually deparaffinize and rehydrate tissue sections mounted on glass slides. For better adhesion of tissue to slide, use of coated or silanized slides is recommended.
5. Immerse room temperature tissue sections into the pre-heated EnVision™ FLEX Target Retrieval Solution (working solution), and incubate for 20-40 minutes at 97 °C. The optimal incubation time should be determined by the user.
6. Leave the sections to cool in PT Link to 65 °C.
7. Remove each slide rack from PT Link tank and immediately dip slides into a jar/tank (e.g., PT Link Rinse Station, Code PT109) with diluted, room temperature EnVision™ FLEX Wash Buffer (20x).
8. Leave slides in the diluted, room temperature EnVision™ FLEX Wash Buffer (20x) for 1-5 minutes.
9. Place slides on a Dako Autostainer/Autostainer Plus instrument and proceed with staining. The sections should not dry out during the treatment or during the immunohistochemical staining procedure.
10. After staining the slides should be mounted using aqueous or permanent mounting medium.

NOTE: When used in PT Link for HIER procedure, the diluted Target Retrieval Solution and diluted Wash Buffer can be used 3 times within a 5 day period if stored at room temperature.

HIER procedure using Coplin jars:

After deparaffinization and hydration to buffer (water), tissue sections should be subjected to heat-induced epitope retrieval (HIER):

1. Place staining jars containing diluted EnVision™ FLEX Target Retrieval Solution (50x) in water bath.
2. Heat water bath and jars filled with EnVision™ FLEX Target Retrieval Solution (working solution) to 95–99 °C. Cover jars with lids to stabilize the temperature and avoid evaporation.
3. Immerse the room temperature sections in the preheated EnVision™ FLEX Target Retrieval Solution (working solution) in the staining jars.
4. Bring temperature of the water bath and EnVision™ FLEX Target Retrieval Solution back to 95–99 °C. Incubate for 20 (±1) minutes at 95–99 °C.
5. Remove the entire jar with slides from the water bath. Allow slides to cool in the EnVision™ FLEX Target Retrieval Solution for 20 (±1) minutes at room temperature.
6. Decant the EnVision™ FLEX Target Retrieval Solution and rinse sections with the diluted, room temperature EnVision™ FLEX Wash Buffer (20x) for 1-5 minutes.
7. Place the rack with rinsed sections in a Dako Autostainer/Autostainer Plus instrument. Immediately begin the pre-programmed staining run.

NOTE: When used in Coplin jars HIER procedure, the diluted Target Retrieval Solution and diluted Wash Buffer can be used only 1 time (single use).

A few epitopes do not tolerate the HIER and some require enzymatic pre-treatment. Please refer to the package insert for the individual Dako primary antibody.

The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended.

B. Staining procedure

Dako FLEX Ready-to-Use Primary Antibodies can be used with Dako EnVision™ FLEX+ detection system, applied on formalin-fixed, paraffin-embedded tissue sections.

A dilution guideline for Dako concentrated primary antibodies is provided in the package insert of the concentrated primary antibody.

The concentrated primary antibodies should be diluted in EnVision™ FLEX Antibody Diluent (DM830, Code K8006).

The staining steps and incubation times are pre-programmed into the software of Dako Autostainer/Autostainer Plus instruments, using the following protocols:

Template protocol: FLEX_200 (200 µL dispense volume) or FLEX_300 (300 µL dispense volume)

Autoprograms for staining runs:

FLEX (FLEX protocol)

FLEX_Mo (FLEX+ Mouse (LINKER) protocol)

FLEX_Rb (FLEX+ Rabbit (LINKER) protocol)

The Auxiliary step should be set to “rinse buffer” in staining runs with ≤10 slides. For staining runs with ≥10 slides the Auxiliary step should be set to “none”. This ascertains comparable wash times.

All incubation steps should be performed at room temperature. For details, please refer to Operator's manual for the dedicated instrument. If the protocols are not available on the used Dako Autostainer/Autostainer Plus instrument, please contact Dako Technical Services

Optimal incubation times of primary antibody and EnVision™ FLEX /HRP are dependent on the applied primary antibody. Please refer to the package insert for the individual Dako primary antibody. The user must verify the applied protocol.

When the staining procedure is completed, the specimens must be mounted. It is recommended to perform dehydration, clearing and permanent mounting when slides have been treated with the 3-in-1 specimen preparation procedure. Slides treated with HIER procedure can be mounted using aqueous or permanent mounting medium. For aqueous mounting, mounting media such as Dako Glycergel™ Mounting Medium, Code C0563, or Faramount Aqueous Mounting Medium, Code S3025, is recommended.

Quality Control

Each staining run should include a known positive control specimen to ascertain a proper performance of all the applied reagents. If the positive control specimen fails to demonstrate positive staining, labeling of test specimens should be considered invalid.

A negative control reagent should be used with each specimen to identify any non-specific staining. If non-specific staining cannot be clearly differentiated from the specific staining, the labeling of the test specimen should be considered invalid.

Interpretation of Results

The diaminobenzidine-containing Substrate Working Solution gives a brown color at the site of the target antigen recognized by the primary antibody. The brown color should be present on the positive control specimen at the expected localization of the target antigen. If non-specific staining is present, this will be recognized as a rather diffuse, brown staining on the slides treated with the negative control reagent. Nuclei will be stained blue by the hematoxylin counterstain.

Limitations

1. For automated use on Dako Autostainer Instruments only. When performing manual stainings, the optimal procedure must be determined and verified by the user, unless specified in the Package Insert for the specific Dako Primary Antibody
2. For use on formalin-fixed, paraffin-embedded tissue sections only. When using frozen sections or cell smears, the optimal procedure must be determined and verified by the user, unless specified in the Package Insert for the specific Dako Primary Antibody. Changes in the procedure may have impact on the number of tests provided with the kit.

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Dako North America, Inc.
Dako Denmark A/S